

Water treatment system



Water that is free from impurities is essential for accurate measurement of acoustic fields and IEC TR62781 provides many recommendations on this matter. Whilst being cheap and readily available, tap water is not ideal for conducting acoustic measurements and should be conditioned before being used for this purpose. The Water treatment system is designed to address these issues. Water that is output from the system will have been subject to:

- Degassing
- Two-stage particulate filtration
- UV filtration

INTRODUCTION

Many ultrasonic measurements require that water is processed to remove water-borne sources of experimental uncertainty. This is discussed at length in IEC TR 62781 ed1.0. Ultrasonics – Conditioning of water for ultrasonic measurements, Geneva Switzerland: International Electrotechnical Commission.

Dissolved gasses can be a significant cause of cavitation and this is the biggest single risk to equipment when measuring high amplitude/high intensity ultrasound fields. Particulates can act as scatters, as food sources for water-borne biological contaminants or as nucleation sites for cavitation. Biological activity within water tank can lead to build up of contaminant on any item placed within the tank and, in extreme cases, are a health risk to the operator. All of these issues are addressed by the Water treatment system.

Many customers already have a source of de-ionised/de-mineralised water on-site. Therefore the Water treatment system from Precision Acoustics Ltd does not include de-ionisation capability. If this is required, customers are recommended to source a supply of de-ionised water locally and alternatives include the use of either ion-exchange resins, reverse-osmosis systems or double/triple distillation of water.

The Water treatment system is designed to be plumbed in line with a water tank system. It is fitted with a timer circuit to allow water to be automatically treated at a time convenient to the user (for example, in the early hours of the morning when the remainder of the ultrasonic test facility is not in use).

TECHNICAL SPECIFICATION

| | |
|-------------------------------|---|
| UV Filtration | Power: 14W Wavelength: 253.7 nm Dosage: >16 mJ/cm ² Lamp lifetime: ≈ 9000 hours |
| Particulate Filtration | Stage 1: 25 micron Stage 2: sub 1micron |
| Degassing | Membrane degassing contactor capable of achieving dissolved O ₂ as low as 2.5 ppm |
| Programmable Timer | 7-day programmable with manual override |
| Maximum operation temperature | 40 °C |
| Water treatment rate | 5-6 litres/minute |

CONSUMABLES

A kit of replacement consumables is available from Precision Acoustics Ltd and includes replacement particulate filters and UV lamp. The rate at which filters become clogged is entirely dependent upon the quality of incoming water. Highly contaminated water may require much more replacement of filter cartridges than cleaner (e.g. distilled) water sources.

All information is based on results gained from experience and tests, and is believed to be accurate but is given without acceptance of liability for loss or damage attributable to reliance thereon as conditions of use lie outside the control of Precision Acoustics Ltd.